

Sequence Processing Software User Manual

Table of Contents

1. Introduction
2. System Requirements and Installation
3. Launching the Application
4. Interface Overview
5. Detailed Module Guide
 - 5.1 Module 1: Conserved Region Analysis
 - 5.2 Module 2: Sequence Export
 - 5.3 Module 3: Difference Calculation
 - 5.4 Module 4: Identity Calculation
 - 5.5 Module 5: Sequence Length Calculation
 - 5.6 Module 6- 1: Fetch Genome from NCBI
 - 5.7 Module 6- 2: Fetch CDS from NCBI
 - 5.8 Module 7: Sequence Optimization (SeqRefine)
 - 5.9 Module 8: Primer Design
 - 5.10 Module 9: 1D/2D Code Generation
6. Common Operations
 - 6.1 Language Switching
 - 6.2 Module Switching
 - 6.3 File and Folder Browsing
7. Output Description
8. Troubleshooting
9. Technical Support and Feedback

1. Introduction

Welcome to the **Sequence Processing Software**, an integrated GUI tool for conserved region analysis, sequence export, difference and identity calculations, length statistics, NCBI data retrieval, sequence optimization, primer design, and 1D/2D code (barcode/QR code) generation. This manual provides step- by- step instructions for

installation, interface navigation, module operation, and troubleshooting to help you get up and running quickly and efficiently.

2. System Requirements and Installation

1. **Operating Systems:** Windows 10/11, macOS 10.15+ , Linux (with PyQt5 support)
2. **Python Version:** Python 3.8 or higher
3. **Required Libraries:**
 - PyQt5
 - Biopython
 - numpy
 - pandas
 - scikit- learn
 - seaborn
 - matplotlib
 - python- docx
 - qrcode
 - pillow
4. **Installation:**

```
pip install PyQt5 biopython numpy pandas scikit-learn seaborn matplotlib python-docx qrcode pillow
```
5. **Program Files:**
 - Place module9.py and 软件图标.png in the same folder. Ensure the icon file is accessible by the application.

3. Launching the Application

1. Open a terminal or command prompt.
2. Change directory:

```
cd /path/to/software/folder
```
3. Run the program:

```
python module9.py
```
4. The main window will appear, ready for input.

4. Interface Overview

- **Top Bar:** Language selector (English/中文) on the left and Module selector (Module 1–9) in the center.
- **Main Panel:** Displays input controls for the selected module using a stacked widget.
- **Run Button:** Each module panel has its own “Run” button at the bottom.
- **Styling:** Gradient background, Arial font, and a layered color scheme for controls.

5. Detailed Module Guide

Each module’s section below describes its purpose, required inputs, usage steps, and outputs.

5.1 Module 1: Conserved Region Analysis

Purpose: Identify highly conserved segments across multiple FASTA sequences.

Inputs:

- **Input Folder:** Directory containing FASTA files (.fasta/.fas).
- **Output Folder:** Directory for saving results.
- **Minimum Length (bp):** Shortest conserved segment to report.
- **Maximum Length (bp):** Longest conserved segment to report.
- **Gap Threshold (0.0–1.0):** Maximum allowed gap proportion per column (default 0.1).
- **Region Threshold (0.0–1.0):** Conservation frequency cutoff per position (default 0.8).

Steps:

1. Click **Browse...** next to “Select folder containing FASTA files” and choose your input folder.
2. Click **Browse...** next to “Select folder to save Conserved Region results” and choose an output folder.
3. Enter numeric values for minimum length, maximum length, gap threshold, and region threshold.
4. Click **Run Module 1**.

Outputs:

- For each FASTA file, a conserved_sequences_<filename>.txt file listing each conserved fragment, its start position, validation flag, and overall SVM cross- validation accuracy.

5.2 Module 2: Sequence Export

Purpose: Split a combined FASTA text file into test and conserved segments for downstream analysis.

Inputs:

- **Input FASTA File:** A .txt file formatted in FASTA style.
- **Output Folder:** Directory for saving split segments.
- **Split Index:** Integer index (e.g. 1) indicating how many sequences are “test” sequences; remaining are conserved.

Steps:

1. Browse and select the input FASTA .txt file.
2. Choose an output folder.
3. Enter the split index.
4. Click **Run Module 2**.

Outputs:

- A set of _segment_<n>.txt files, where each file contains all sequences for that segment.

5.3 Module 3: Difference Calculation

Purpose: Compute average nucleotide- level difference between test sequences and species segments, output distribution tables and heatmaps.

Inputs:

- **Test Sequence File:** The .txt file with test sequences.
- **Species Folder:** Folder containing _segment_<n>.txt files.
- **Output Excel File:** Path for saving results (.xlsx).

Steps:

1. Select the test sequence .txt file.
2. Select the species segment folder.
3. Specify an output .xlsx file path.
4. Click **Run Module 3**.

Outputs:

1. **Excel Workbook** with three sheets:
 - Avg Nuc Diff (%): Average difference per segment.
 - Species Diff (%) (Prop): Proportion distribution for thresholds 1–99%.
 - Species Diff (%) (Num): Count distribution for bases thresholds.

2. **Heatmaps:** Saved as <prefix>_species_diff_prop_heatmap.svg/png and <prefix>_species_diff_num_heatmap.svg/png.

5.4 Module 4: Identity Calculation

Purpose: Compute average nucleotide- level identity against species consensus, with distribution details and heatmaps.

Inputs: Same as Module 3.

Steps:

- Provide test sequence and species folder paths.
- Enter output Excel path.
- Click **Run Module 4**.

Outputs:

1. **Excel Workbook** with:
 - Avg Nuc Identity (%)
 - Species Identity (%) (Prop)
 - Species Identity (%) (Num)
2. **Heatmaps:** <prefix>_species_identity_prop_heatmap.svg/png, <prefix>_species_identity_num_heatmap.svg/png.

5.5 Module 5: Sequence Length Calculation

Purpose: Calculate the length of each sequence in a FASTA file.

Inputs:

- **Input FASTA File:** .fasta or .fas.
- **Output Excel File:** Path to .xlsx.

Steps:

1. Browse for the FASTA file.
2. Specify output Excel path.
3. Click **Run Module 5**.

Outputs:

1. An Excel sheet with columns: Title and Length.

5.6 Module 6-1: Fetch Genome from NCBI

Purpose: Batch- download complete genome FASTA for given accession numbers via Entrez.

Inputs:

- **Accession List File:** .txt file listing one accession number per line.
- **Output Folder:** Destination for .fasta files.

Steps:

1. Select the accession list file.
2. Choose an output folder.
3. Click **Run Module 6-1**.

Outputs:

- Individual <accession>.fasta files for each number in the list.

5.7 Module 6-2: Fetch CDS from NCBI

Purpose: Batch- download GenBank records, extract CDS features, and save them as text.

Inputs: Same as Module 6-1.

Steps:

1. Provide the accession list file.
2. Select an output directory.
3. Click **Run Module 6-2**.

Outputs:

- <accession>_cds.txt files containing CDS feature entries.

5.8 Module 7: Sequence Optimization (SeqRefine)

Purpose: Remove degenerate bases and optionally rename sequence titles in FASTA files.

Inputs:

1. **Input Folder:** Contains .fasta or .fas files.
2. **Modify Titles:** Check to use a custom title.
3. **New Title:** Text (without >) if modifying titles.
4. **Rename to Filename:** Check to rename each sequence header to its file basename.

Steps:

- Choose the FASTA folder.
- (Optional) Tick **Modify Titles** and enter a new title.

- (Optional) Tick **Rename to Filename**.
- Click **Run Module 7**.

Outputs:

- FASTA files overwritten in place: degenerate bases removed and headers updated accordingly.

5.9 Module 8: Primer Design

Purpose: Automatically design forward and reverse primers, calculate GC content and melting temperature.

Inputs:

- **Input FASTA File:** .fasta or .fas.
- **Output Excel File:** Path for .xlsx.
- **Primer Length:** Integer (default 18 bp).

Steps:

- Browse to select the input FASTA.
- Specify output Excel path.
- Enter desired primer length.
- Click **Run Module 8**.

Outputs:

1. Excel table listing region names, forward/reverse primer sequences, GC%, and Tm.

5.10 Module 9: 1D/2D Code Generation

Purpose: Convert text or DOCX content to stylized barcodes/QR codes and generate a Word document with base- replacements.

Inputs:

- **Input File:** .txt or .docx.
- **Output File:** .docx for final document.

Steps:

1. Select the input text or Word file.
2. Set the output .docx path.
3. Click **Run Module 9**.

Outputs:

- Stylized PNG barcode/QR images saved to the output folder.
- A Word document containing replaced base symbols and embedded codes.

6. Common Operations

6.1 Language Switching

Use the top- left combo box to switch between **English** and 中文—all UI text updates instantly.

6.2 Module Switching

Select any module (1–9) from the top- center combo box; the main panel will display that module’s inputs.

6.3 File and Folder Browsing

Click each panel’s **Browse...** button to open file/folder dialogs. Selected paths populate the adjacent text fields.

7. Output Description

- All generated files (Excel, TXT, PNG, SVG, DOCX) appear in user- specified directories.
- Filenames include module- specific prefixes/suffixes.
- Heatmaps: files ending with `_prop_heatmap` (percentage) or `_num_heatmap` (count).

8. Troubleshooting

1. **Missing Dependencies:** Ensure all Python packages are installed.
2. **Invalid Paths:** Verify file/folder existence and read/write permissions.
3. **Excel Export Errors:** Close open Excel instances or change the output filename.
4. **NCBI Fetch Failures:** Check Internet connectivity and set `Entrez.email` in the code.
5. **GUI Unresponsive:** Restart the application; confirm Python and PyQt5 compatibility.

9. Technical Support and Feedback

For questions or feature requests, contact our support team:

- **Email:** gs2022@hnu.edu.cn

Thank you for using the Sequence Processing Software. We wish you success in your research!